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Drug Class	Drug	Monitoring/Measuring Test	Optimal or Desired Time of Blood Collection	Notes
Heparin	Unfractionated	APTT Anti-Xa	6 hours after dose initiation or adjustment	<ul style="list-style-type: none"> • Anti-Xa measurements also acceptable. • Circuit anticoagulation may require more frequent monitoring. • Circuit anticoagulation may require higher dosing that cannot be measured by APTT, so ACT may be the optimal test. • Hybrid or LMWH calibrated anti-Xa are acceptable. • Pentasaccharide-calibrated anti-Xa reported in mg/L or mg/dL.
	Low-molecular-weight	Anti-Xa	4 hours after third dose	
	Pentasaccharide	Anti-Xa	3 hours after dose	
Vitamin K Antagonists	Antagonist—oral	PT/INR (baseline PT/INR should be collected prior to initiation of therapy)	First INR is with 12–24 hours of first dose	<ul style="list-style-type: none"> • If outside target, consider VKORC mutation. • If within target, dose is probably acceptable. • If lower than target, probably need to adjust dose higher. • Patient may have additional reversal agents. • More-frequent monitoring may be necessary if patient is bleeding.
	(reversal therapy)	PT/INR	IV: 12 hours Oral: 12–24 hours	
Antithrombotics	Antiplatelet—oral	Platelet-function tests	1 week after initiation of therapy	<ul style="list-style-type: none"> • No consensus on method for measuring or whether necessary. • Methods include PFA System, ULTEGRA, and traditional platelet aggregation. • POC or platelet-function methods, including ULTEGRA, TEG, or ROTEM-based methods. Traditional platelet aggregation studies may also be suitable. • No consensus or recommendations for monitoring available. • Platelet-function tests (if concomitant therapy, then ACT) • 30 minutes for checking catheter function. • No consensus or recommendations for monitoring available.
	Antiplatelet— intravenous	Platelet-function tests (if concomitant therapy, then ACT)	Drug-dependent Within 0.5–2 hours	
	Fibrinolytics	TT, FBG, XDP, FDP	10–15 minutes after completion of infusion	
	Anti(fibrin)olytics	PT/INR, APTT, ACT	Max. concentration at 2 hours after initiation of infusion	
	DTI	APTT, anti-II, ECA, ECT (may require ACT if high dose for PCI), dilute thrombin time	HIT treatment: 4–6 hours PCI infusion: 5 and 45 minutes	<ul style="list-style-type: none"> • Drug labeling indicates APTT or ACT, but other more-specific testing may be desirable, especially when therapeutic target is not achieved.
DOAC	IIa	ECT, ECA, IIa, dTT	Trough samples (5–30 minutes prior to next dose)	<ul style="list-style-type: none"> • Collect just before next dose. • If peak samples are desired, then usually 2–3 hours after dose. • Collect just before next dose. • If peak samples are desired, then usually 2–3 hours after dose. • Specific drug-calibrated anti-Xa. • No consensus or recommendations for monitoring. • No consensus or recommendations for monitoring. • The drug is continuously infused. Once infusion has stopped, DOAC levels may rise. • Specific DOAC-calibrated anti-Xa. • Note that newer replacement therapies for hemophilia using modified (PEGylated, albumin-fused, etc) factor replacement may require special methods. • Single-stage clotting assay or chromogenic assay may be preferred. • Drug delivery is either by nasal spray or infusion. • No consensus or recommendations for monitoring. • For APCC, no relationship to decreasing clotting times and bleeding outcomes. • No consensus or recommendations for monitoring. • FBG, FVIII, FXIII, or other factors may be assessed if replacing for that purpose.
	Xa	Drug-calibrated anti-Xa	Trough samples (5–30 minutes prior to next dose)	
	Reversal—dabigatran	ECT, ECA, anti-IIa	10–15 minutes after completion of infusion	
	Reversal—rivaroxaban/apixaban	Anti-Xa	4 hours after infusion to reassess anti-Xa	
Factor Replacement	Replacement therapy (hemophilia A or B)	Factor level	Physician-guided. PK studies may be required. Consider baseline, 30 min, 60 min, 2, 4, 8, 12, and 24 hours for PK time periods.	<ul style="list-style-type: none"> • Note that newer replacement therapies for hemophilia using modified (PEGylated, albumin-fused, etc) factor replacement may require special methods. • Single-stage clotting assay or chromogenic assay may be preferred. • Drug delivery is either by nasal spray or infusion. • No consensus or recommendations for monitoring. • For APCC, no relationship to decreasing clotting times and bleeding outcomes. • No consensus or recommendations for monitoring. • FBG, FVIII, FXIII, or other factors may be assessed if replacing for that purpose.
	DDAVP/vasopressin	VWF, FVIII	Baseline, 30 minutes, 2, 4, and 6 hours post-drug delivery	
	PCC or APCC	PT/INR, APTT	Baseline (pretreatment) and 15–30 minutes after completion of administration	
	rVlla	PT/INR	10–15 minutes after infusion	
	FFP, cryoprecipitate, or other plasma-based products	PT/INR, APTT	Baseline (preinfusion) and 30 minutes	

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
Odběr vzorku

Order of Draw	Vacuum Tube Type	Color of Cap	Test Type
1	Blood culture tubes	Color varies	Blood culture
2	Sodium citrate (3.2%)	Light blue	Coagulation
3	Glass (no activator)	Red	Chemistry, immunoassays
4	SST	Gold or red/black	Chemistry, serology
5	Trace elements (no preservative)	Royal blue	Trace element toxicology
6	Sodium or lithium heparin	Green	Chemistry
7	EDTA	Lavender or pink	Hematology, blood bank
8	Sodium fluoride	Gray	Glucose
9	ACD	Yellow	Blood bank HLA testing
10	QUANTIFERON-TB	Gold	TB testing

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Coagulation activation² (partial clotting)

Assay	Effect
PT sec, PT INR	↑
PT %	↓
APTT sec	↑ ↓
TT sec	↑
FII, FV, FVIII %	↓
FVII %	↑
D-Dimer mg/l	↑
VWF %	↓

Sample dilution (excess CaCl₂-Citrate binding)⁴

Assay	Effect
PT sec, PT INR	↑
PT %	↓
APTT sec	↑
Factors %	↓
TT sec	↑
D-Dimer mg/l	↓

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Pre-analytical Variables In Routine Coagulation Testing: Setting the Stage for Accurate Results

Robert C. Gosselin^a and Richard A. Marlar^b

^aHemophilia Treatment Center, University of California, Davis Health System, Sacramento, CA, U.S.

^bDepartment of Pathology, University of New Mexico Health Sciences Center Albuquerque, NM, U.S.

Introduction

Many pre-analytical variables (PAVs) may affect the results of routine coagulation assays. To improve the precision and accuracy of laboratory testing, it is critical to identify these variables and realize their potential impact.^{1,2} Additionally, advances in laboratory instrumentation have improved the reproducibility and sensitivity of the analytical phase, therefore creating greater dependence on specimen integrity.^{3,5} The activated partial thromboplastin time (APTT) and prothrombin time (PT) determinations are among the most frequently ordered screening tests in the clinical laboratory. These assays are used in the evaluation of a wide variety of clinical conditions, for either diagnostic or monitoring purposes. These screening tests also form the basis of many special coagulation tests, such as factor assays and proteins C and S activity assays.

Since the introduction of coagulation assays, efforts have been made not only to automate these assays, but to better standardize testing, thus providing more-accurate results to aid in clinical assessment.⁴⁻⁶ PAVs pertaining to routine coagulation testing can be classified into three major categories:

1. Specimen collection (including patient selection)
2. Specimen transportation and stability
3. Specimen processing and storage

Within each of these categories, there are a number of individual variables, each of which may have a major impact on testing. There are also a number of variables related to the analysis of the specimen, many of which are dependent on the reagent(s) and instrumentation, but variables associated with analysis are beyond the scope of this document.

Many standards for testing in the general clinical laboratory and specifically in the coagulation laboratory have been developed in an effort to improve precision and accuracy.⁷⁻¹⁰ The Clinical and Laboratory Standards Institute (CLSI) is the primary organization for clinical laboratory standards and guidance documents in the

United States, although their guidelines are also referenced and used internationally.⁴ The CLSI was established in 1968 as a group of individuals representing industry, government, and professionals dedicated to the development of standards and guidance documents for clinical laboratory testing. Despite CLSI's existence, there is still a lack of practice standardization among clinical laboratories regarding specimen collection, storage, and processing for coagulation testing. Some of the procedures in practice today are apparently founded on tradition, while others are based on CLSI guidelines, with and without significant published or supporting evidence. As a result, a number of problems, inconsistencies, and erroneous results can still arise based on pre-analytical processing of the specimen, and these discrepancies may be associated with disastrous outcomes. This manuscript will review the pre-analytical variables and some analytical variables of laboratory-based coagulation testing. Point-of-care (POC) devices, which may use native or anticoagulated whole blood, are also affected by PAVs but will have limited focus in this document. Recommendations for the proper methods based on published and new data will be presented, together with recommendations for converting to these methods.

Pre-analytical Variables In Routine Coagulation Testing: Setting the Stage for Accurate Results

Gosselin RC, University of California, Davis Health System, Sacramento, California, U.S.

Marlar RA, University of New Mexico Health Sciences Center Albuquerque, New Mexico, U.S.